

## **ENZYMATIC TEST KIT FOR THE DETERMINATION OF CITRIC ACID IN GRAPE JUICE AND WINE**

### **PRODUCT**

Product no. 4A126, **for 30 tests, for *in vitro* use only**

### **PRINCIPLE OF MEASUREMENT**

Citric acid may be used at the final stages of winemaking to make minor adjustments to acid levels without affecting the bi-tartrate stability of the wine. It is determined enzymatically according to the following equations:



In the presence of enzymes MDH and LDH, both the oxaloacetate and its decarboxylation product pyruvate, are reduced by NADH to malate and lactate respectively.



The amount of NADH oxidised is measured at 340 nm and is stoichiometrically related to the amount of citrate present.

### **CONTENTS**

The kit includes the following reagents:

| Reagent No. | Reagent  | Preparation   | Quantity   | Stability  |
|-------------|----------|---|------------|--|
| 1           | Buffer   | Nil   | 33 mL      | Reagents are stable for 12 months at 4°C or until the kit's expiry date, whichever occurs first. Reagent 2 and 4 are stable for 1 month (Reagent 2) and 2 months (Reagent 4) at 4°C once dissolved or until the kit's expiry date, whichever occurs first. |
| 2           | NADH     | Add 1.7 mL of distilled water to either bottle as required, mix to dissolve | 2 x 1.7 mL |  |
| 3           | MDH/LDH  | Mix gently by inversion before use  | 0.7 mL     |  |
| 4           | CL       | Add 0.35mL of distilled water to either bottle as required, mix to dissolve | 2 x 0.35mL |  |
| 5           | Standard | Nil   | 3.3 mL     |  |

The shelf life of Reagent 1 can be extended by placing aliquots in a freezer.

Do not freeze reagents 2, 3 or 4.

Failure to store reagents at the recommended temperature will reduce their shelf life.

For concentration of Standard, refer to label on bottle.

### **SAFETY**

- Wear safety glasses
- Do not ingest Buffer or Standard as they contain sodium azide as a stabilizer

### **PROCEDURE**

Operating Parameters

|                         |  |
|-------------------------|--|
| Wavelength              | 340 nm   |
| Cuvettes                | 1cm, quartz, silica, methacrylate or polystyrene |
| Temperature             | 20 – 25°C  |
| Final volume in cuvette | 3.14 mL  |
| Zero                    | against air without cuvette in light path        |

### SAMPLE PREPARATION

Samples should be diluted to ensure concentration in the assay solution is no more than 0.5 g/L. For most samples, a 1 in 2 dilution with distilled water should be sufficient. For samples containing between 1 g/L to 2.5 g/L of citric acid, a 1 in 5 dilution would be appropriate. Ideally,  $A_1$  should lie between 0.90 – 1.20 absorbance units.

Red wines or highly coloured undiluted juice samples require decolourisation. To decolourise, add approximately 0.1 g of PVPP to 5 mL of sample in a test tube. Shake well for about 1 minute. Clarification is achieved by settling or filtering through Whatman No. 1 filter paper.

### SAMPLE ANALYSIS

a. Pipette the following volumes of reagents into the cuvettes:

| Reagent         | Blank                  | Standard               | Sample                 |
|-----------------|------------------------|------------------------|------------------------|
| 1. Buffer       | 1.00 mL (1000 $\mu$ L) | 1.00 mL (1000 $\mu$ L) | 1.00 mL (1000 $\mu$ L) |
| 2. NADH         | 0.10 mL (100 $\mu$ L)  | 0.10 mL (100 $\mu$ L)  | 0.10 mL (100 $\mu$ L)  |
| Distilled water | 2.00 mL (2000 $\mu$ L) | 1.80 mL (1800 $\mu$ L) | 1.80 mL (1800 $\mu$ L) |
| Sample/Standard |                        | 0.20 mL (200 $\mu$ L)  | 0.20 mL (200 $\mu$ L)  |
| 3. MDH/LDH      | 0.02 mL (20 $\mu$ L)   | 0.02 mL (20 $\mu$ L)   | 0.02 mL (20 $\mu$ L)   |

b. Mix well by gentle inversion and read absorbances,  $A_1$ , after 5 minutes.

c. Pipette the following reagent into the cuvettes:

|       |                      |                      |                      |
|-------|----------------------|----------------------|----------------------|
| 4. CL | 0.02 mL (20 $\mu$ L) | 0.02 mL (20 $\mu$ L) | 0.02 mL (20 $\mu$ L) |
|-------|----------------------|----------------------|----------------------|

d. Mix well by gentle inversion and read absorbances,  $A_2$ , after 25 minutes.

### CALCULATIONS\*

1. Calculate the Net Absorbance for the Blank, Sample and Standard:

$$\text{Net Absorbance, } A_N = A_1 - A_2$$

2. Calculate the Corrected Absorbance by subtracting the Net Absorbance for the Blank from the Net Absorbance for the Sample.

$$\text{Sample Corrected Absorbance, } A_C = \text{Sample } A_N - \text{Blank } A_N$$

3. Do the same for the Standard by substituting the Standard absorbance values in place of the Sample absorbance values.

4. Calculate the Citric acid concentration as follows;

$$\text{Citric acid (g/L)} = A_C \times 0.4787 \times \text{Dilution Factor}$$

\*A calculation spreadsheet is available for download at the following locations:

Australia based users

<https://winechek.com/calculation-worksheets/>

Users outside of Australia

<http://www.vintessential.com.au/resources/calculation-worksheets/>

### REFERENCES

1. OIV, 2018, Compendium of international methods of wine and must analysis. *International Organisation of Vine and Wine*, Vol 1: Paris, France, pp. OIV-MA-AS313-09.